


PROCESS FILTER HOUSINGS & EQUIPMENTS




ON LINE HOUSINGS

ON LINE HOUSINGS are designed for vertical online mounting primarily for air vents on the tanks. These filter housings have good flow characteristics & can be rinsed easily. These type of filter housings are mostly used for online filtration of products in Pharmaceuticals, Food & Beverages and Biotechnology industries.



MULTIROUND ONLINE FILTER HOUSING



SINGLE ROUND ONLINE FILTER HOUSING

| PART NO. | "A" (IN MM) | | | | "B" (IN MM) | "C" (IN MM) ID | CONNECTION | | | | BOX SIZE IN INCH (LxWxH) | | | |
|---------------|-------------|-----|-----|------|----------------|----------------------|-----------------|------------|-------------|----------|-------------------------------|----------|----------|--|
| | 10" | 20" | 30" | 40" | | | INLET OUTLET | VENT TC | DRAIN TC | 10" | 20" | 30" | 40" | |
| | | | | | | | | | | | | | | |
| A0401-01XXXBC | 273 | 528 | 795 | 1048 | 75 | 98 | 1.5" | 1/2" | 1/2" | 06x06x20 | 06x06x30 | 06x06x40 | 06x06x50 | |
| A0401-04XXXBS | 300 | 550 | 800 | 1050 | 330 | 214 | 1.5" | 1/2" | 1/2" | 20x20x34 | 20x20x44 | 20x20x54 | 20x20x64 | |
| A0401-07XXXBS | 300 | 550 | 800 | 1050 | 340 | 250 | 2" | 1/2" | 1/2" | 20x20x35 | 20x20x45 | 20x20x55 | 20x20x65 | |
| A0401-09XXXBS | 300 | 550 | 800 | 1050 | 350 | 300 | 2" | 1/2" | 1/2" | 22x22x36 | 22x22x46 | 22x22x56 | 22x22x66 | |
| A0401-12XXXBS | 300 | 550 | 800 | 1050 | 365 | 350 | 2" | 1/2" | 1/2" | 24x24x36 | 24x24x46 | 24x24x56 | 24x24x66 | |
| A0401-16XXXBS | 300 | 550 | 800 | 1050 | 400 | 400 | 2" | 1/2" | 1/2" | 26x26x37 | 26x26x47 | 26x26x57 | 26x26x67 | |
| A0401-21XXXBS | 300 | 550 | 800 | 1050 | 405 | 450 | 2" | 1/2" | 1/2" | 28x28x37 | 28x28x47 | 28x28x57 | 28x28x67 | |
| A0401-25XXXBS | 300 | 550 | 800 | 1050 | 440 | 500 | 4" | 1/2" | 1" | 30x30x39 | 30x30x49 | 30x30x59 | 30x30x69 | |
| A0401-31XXXBS | 300 | 550 | 800 | 1050 | 450 | 550 | 4" | 1/2" | 1" | 32x32x40 | 32x32x50 | 32x32x60 | 32x32x70 | |
| A0401-37XXXBS | 300 | 550 | 800 | 1050 | 460 | 600 | 4" | 1/2" | 1" | 34x34x40 | 34x34x50 | 34x34x60 | 34x34x70 | |
| A0401-40XXXBS | 300 | 550 | 800 | 1050 | 470 | 650 | 6" | 1/2" | 1.5" | 36x36x41 | 36x36x51 | 36x36x61 | 36x36x71 | |
| A0401-47XXXBS | 300 | 550 | 800 | 1050 | 480 | 700 | 6" | 1/2" | 1.5" | 38x38x41 | 38x38x51 | 38x38x61 | 38x38x71 | |
| A0401-54XXXBS | 300 | 550 | 800 | 1050 | 490 | 750 | 6" | 1/2" | 1.5" | 40x40x42 | 40x40x52 | 40x40x62 | 40x40x72 | |

XXX :- Length of Cartridge (10",20",30",40") & End Connections

ORDERING GUIDE:

