


PROCESS FILTER HOUSINGS & EQUIPMENTS




T-LINE HOUSINGS

T-LINE HOUSINGS are designed for providing a robust solution for critical filtration applications. These housings are often more convenient as it is free standing & the cartridge can be changed without breaking inlet & outlet connections. T-line filter housings are used for sterile filtration, separation of micro-particles and many more applications. Chemical, Pharmaceutical, Biotechnology industries use these type of filter housings.



MULTIROUND T LINE FILTER HOUSING



SINGLE ROUND T LINE FILTER HOUSING

PART NO.	"A" (IN MM)				"B"	"C"	"D" (IN MM)	CONNECTION								BOX SIZE IN INCH (LxWxH)			
	10"	20"	30"	40"				INLET OUTLET	PRESSURE GAUGE		VENT		DRAIN/SAMPLING		10"	20"	30"	40"	
									TC	BSP	TC	BSP	TC	BSP					
A0402-01XXXBC	273	528	795	1048	250	250	98	1"	1.5"	1/4"	1/2"	1/4"	1/2"	1/4"	15x15x29	15x15x39	15x15x49	15x15x59	
A0402-04XXXBS	300	550	800	1050	400	325	214	1.5"	1.5"	1/4"	1/2"	1/4"	1/2"	1/4"	20x20x35	20x20x45	20x20x55	20x20x65	
A0402-06XXXBS	300	550	800	1050	400	340	250	2"	1.5"	1/4"	1/2"	1/4"	1/2"	1/4"	20x20x36	20x20x46	20x20x56	20x20x66	
A0402-08XXXBS	300	550	800	1050	450	370	300	2"	1.5"	1/4"	1/2"	1/4"	1/2"	1/4"	23x23x37	23x23x47	23x23x57	23x23x67	
A0402-10XXXBS	300	550	800	1050	500	380	350	2"	1.5"	1/2"	1/2"	1/2"	1/2"	1/2"	25x25x38	25x25x48	25x25x58	25x25x68	
A0402-16XXXBS	300	550	800	1050	550	390	400	2"	1.5"	1/2"	1/2"	1/2"	1/2"	1/2"	27x27x40	27x27x50	27x27x60	27x27x70	
A0402-20XXXBS	300	550	800	1050	600	395	450	2"	1.5"	1/2"	1/2"	1/2"	1/2"	1/2"	29x29x40	29x29x50	29x29x60	29x29x70	
A0402-24XXXBS	300	550	800	1050	650	450	500	3"	1.5"	1/2"	1/2"	1/2"	1"	1"	31x31x44	31x31x54	31x31x64	31x31x74	
A0402-27XXXBS	300	550	800	1050	700	460	550	3"	1.5"	1/2"	1/2"	1/2"	1"	1"	32x32x46	32x32x56	32x32x66	32x32x76	
A0402-36XXXBS	300	550	800	1050	750	480	600	3"	1.5"	1/2"	1/2"	1/2"	1"	1"	35x35x48	35x35x58	35x35x68	35x35x78	
A0402-40XXXBS	300	550	800	1050	800	500	650	4"	1.5"	1/2"	1/2"	1/2"	1.5"	1.5"	37x37x49	37x37x59	37x37x69	37x37x79	
A0402-45XXXBS	300	550	800	1050	850	520	700	4"	1.5"	1/2"	1/2"	1/2"	1.5"	1.5"	38x38x50	38x38x60	38x38x70	38x38x80	
A0402-50XXXBS	300	550	800	1050	900	550	750	4"	1.5"	1/2"	1/2"	1/2"	1.5"	1.5"	41x41x51	41x41x61	41x41x71	41x41x81	

XXX :- Length of Cartridge (10",20",30",40") & End Connections

ORDERING GUIDE:

