

# FILTRATION IN COATING CHEMICALS



We at MICROFILT INDIA PVT LTD take pleasure to inform you that our company has completed more than two decades in Filtration and Separation Technology catering to critical filtration need in manufacturing of Pharmaceutical, Biotech, coating Chemicals and Food Beverage.



## FEATURE

1. Consistent quality filtration.
2. High temperature filtration
3. High viscosity filtration.

When it comes to filtration of Coating chemical solution we understands the stress and focus that you place the quality of the product that you produce.

Latest technologies permit the possibility of allowing the reduction of thickness of paint layers to very low levels. This put challenge to filtration operation.

## APPLICATION

- Paint
- Resins
- Ink
- Varnish
- Lacquers
- Polymers



## STRATEGY ADOPTED

We have excellent product and process to optimize your process of filtration to remove undesired element. MICROFILT provide MAS (Microfilt automatic self cleaning strainers) filter for trouble free operation. Many time BAG filter is used in combination with MAS. filter is used in combination with MAS.



## MAS

1. Eliminate Use of disposable media
2. Choice of metal cleanable media
3. Mechanically clean in place filters
4. Required lower floor space
5. Continuous and consistence flow filtration

## SUCCESSFUL INSTALLATION

- Microink pvt.ltd.
- Zydex pvt.ltd
- Arkema pvt.ltd.
- Astra chemtech pvt.ltd.
- Grauer &weil pvt.ltd.
- KLJ polymers and plasticizers pvt. Ltd.
- And many more



  
**Let's Colour**



**MICROFILT INDIA PVT LTD.**

**247/2,GIDC Phase II, Umbergaon**

**Dist. Valsad Gujarat-396171**

**(0260)2561392 EXT:34**

**E mail:-marketing@microfiltindia.com**